

stryker®

Neuro Spine ENT

ESSx™ Microdebrider

Your Specialty Is Our Priority.



**Introducing the ergonomically
advanced ESSx shaver.**

ESSx™


Microdebrider

stryker®

Neuro Spine ENT


ESSx Microdebrider

5290-601-200



- **Max RPM in oscillation mode (standard):** 5,000 RPM
- **Straight through suction:** Helps ensure less clogging
- **Improved patented flush feature:** Surgeon can manually unplug hand piece using its own irrigation out flow allowing for quicker removal of clots and decreased OR time
- **Integrated scope cleaning:** Cleans the endoscope allowing for less manual cleaning during the case
- All straight blades are irrigated through the entire channel including the tip
- Preserves mucous membranes
- Improved ergonomics and balance
- Increased cutting efficiency

Precision Cutters



SharpEdge

Part #	Dia.	Length
5290-340-000	4.0 mm	8cm



Aggressive

Part #	Dia.	Length
5290-648-000	4.0 mm	8cm
5290-638-000	3.5 mm	8cm
5290-628-000	2.5 mm	8cm
5290-011-000	4.0 mm	11cm



Aggressive Serrated

Part #	Dia.	Length
5290-645-000	4.0 mm	8cm
5290-635-000	3.5 mm	8cm
5290-638-000	3.5 mm	8cm
5290-628-000	2.5 mm	8cm
5290-011-000	4.0 mm	11cm



Jaguar

Part #	Dia.	Length
5290-540-000	4.0 mm	8cm
5290-530-000	3.5 mm	8cm



Angled Aggressive 60°

Part #	Dia.	Length
5290-760-100	4.0 mm	11cm



Angled Aggressive 40°

Part #	Dia.	Length
5290-740-100	4.0 mm	11cm



Angled Aggressive 12°

Part #	Dia.	Length
5290-744-100	4.0 mm	11cm



Angled Aggressive 40° Convex

Part #	Dia.	Length
5290-740-200	4.0 mm	11cm



Angled Aggressive Serrated 40°

Part #	Dia.	Length
5290-440-100	mm	cm



Angled Aggressive Serrated 60°

Part #	Dia.	Length
5290-640-100	mm	cm

Burs



Hooded Round Bur

Part #	Dia.	Length
5290-647-000	3.0 mm	8cm



Unhooded Round Bur

Part #	Dia.	Length
5290-220-000	2.0 mm	8cm

Accessories



iSwitch Receiver
0277-200-100



CORE Console
5400-050-000



iSwitch Footpedal
0277-100-100



Bi-Directional Footpedal
5100-008-000

Irrigation Cassettes
5290-075-000



Mini-Trephination Set
0290-222-000

The information presented in this brochure is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Products referenced with ™ designation are trademarks of Stryker. Products referenced with ® designation are registered trademarks of Stryker.

Literature Number: 9100-001-091 Rev. A UnDe/P.S.

Copyright © 2007 Stryker
Printed in USA

Stryker® Neuro Spine ENT
4100 East Millham Avenue
Kalamazoo, MI 49001 USA
t: 269 323 7700 f: 269 353 3026
toll free: 800 253 3210

www.stryker.com

