

OtoPore®

# Bioresorbable Ear Dressings



**For outer and middle ear surgery**

Made of bioresorbable synthetic polymer, OtoPore ear dressings provide post-op structural support in various shapes and sizes to meet the needs of different anatomical placements and user preferences. Its body-friendly formula self degrades into natural components and clears through the Eustachian tube or mucosal flow.



## Predictable performance

OtoPore bioresorbable synthetic polyurethane foam is available in different densities, shapes and sizes. It aids the post-op structural support.

- Provides hemostasis through compression
- Minimizes adhesions and edema
- Absorbs blood/fluid without swelling in size
- Separates and supports vital tissue structures
- Clears naturally; atraumatic and comfortable for patients
- Easy to cut and manipulate during placement
- No mixing, moistening or other prep needed

	Ear Cylinder Standard	Ear Cylinder Firm	Ear Square Firm
Order Number (QTY/8 per box)	5400-010-000 (2 x 1cm)	5400-020-000 (2 x 1cm)	5400-020-100 (2cm x 4cm x 1.5mm)
Typical Placement (more than one unit may be applied if needed)	Outer ear	Outer ear	Middle ear
Degradation Time*	7-15 days	7-15 days	14-28 days
Support Time	3-5 days	4-6 days	4-10 days
Compression	36-48 hours	36-48 hours	36-48 hours
Hydrophilic (does not swell)	✓	✓	✓
Material	100% synthetic proprietary polymer	100% synthetic proprietary polymer	100% synthetic proprietary polymer

\*May vary due to amount used, degree of saturation and site

## Ordering is easy


All variations of our bioresorbable solutions are conveniently available through a single source. To order, call Customer Service at **800-253-3210** or contact your Stryker NSE sales representative. For more details visit **[www.stryker.com/bioresorbables](http://www.stryker.com/bioresorbables)**



The information presented in this brochure is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Products referenced with TM designation are trademarks of Stryker. Products referenced with ® designation are registered trademarks of Stryker.

Literature Number: 9100-003-282 Rev None G55

 See OtoPore Instructions for Use (IFU) for complete details.

Copyright © 2015 Stryker  
Printed in USA

Stryker Neuro Spine ENT  
4100 East Mitham Avenue  
Kalamazoo, MI 49001 USA  
toll free: 800 253 3210

[www.stryker.com/NSE](http://www.stryker.com/NSE)